


## HPLC DATA SHEET

1

## アシクロビル (日本薬局方記載条件) Aciclovir (The Japanese Pharmacopoeia)


H140429A

(A) Related substances: Standard solution\*1




	System suitability requirement	Result
Relative standard deviation of the peak area	≤2.0%	1.87%
(Guanine)	<b>≟</b> 2.0 /6	1.07 /0
Relative standard deviation of the peak area (Aciclovir)	<b>≦</b> 1.0%	0.10%
Resolution (1,2)	≧17	42.3

(B) Assay: Standard solution\*1 (Aciclovir 1.0 mg/mL)



HN N Aciclovir

(C) Assay: System suitability solution\*1 (Aciclovir 1.0 mg/mL, Guanine 0.005 mg/mL)



Column : Meteoric Core C18 BIO (2.7 µm, 16 nm)

100 X 4.6 mml.D.

Eluent : phosphate buffer containg 1-decansulfonic acid sodium salt\*2/acetonitrile (1000/40)

<sup>2</sup>Dissolve 1.0 g of containg 1-decansulfonic acid sodium salt and 6.0 g of NaH<sub>2</sub>PO<sub>4</sub> • 2H<sub>2</sub>O in

1000 mL of water, adjust pH 3.0 with H<sub>3</sub>PO<sub>4</sub>

Flow rate : 1.4 mL/min (adjust the flow rate so that the retention time of aciclovir is about 3 min)

Temperature : 20°C

Detection : UV at 254 nm

Injection : 10 µL

Pressure : 26.5 MPa (3840 psi)

(The Japanese Pharmacopoeia 16 th; Related substances, Assay)

\*1 All solutions were prepared from aciclovir and guanine supplied as a reagent for laboratory use.